




# SKYSCAN 1272

- High-Resolution X-ray Microtomograph

## Polymers & Composites


- Resolve fine structures with voxel sizes down to 0.35  $\mu\text{m}$
- Assess microstructural composition, architecture, and porosity
- Quantify defects, fiber thickness, and orientation



Carbon fiber reinforced polymer, after subjection to impact test; 8  $\mu\text{m}$  voxel size, 70 kV, 0.5 mm Al + 0.38 mm Cu filter, 4,720 x 1,088 x 1,200 px – 40 mm sample

## Geology, Oil & Gas


- Measure pore network properties, grain size, and shape
- Calculate distribution of mineral phases in 3D
- Digitize a 3D volume of precious samples, e.g. archeological finds
- Analyze dynamic processes



Sandstone microplug, pore network color-coded according to size  
1  $\mu\text{m}$  voxel size, 80 kV, 0.5 mm Al filter, 1,332 x 2,000 x 1,050 px – 2 mm sample size

## Biomaterials


- Apply in-situ analysis of mechanical properties and dynamic processes
- Quantify porosity, pore network properties, and local thicknesses in 3D
- Detect inhomogeneities and deviations in the printing process



PLA scaffold, central part color-coded for local thickness of struts  
2.5  $\mu\text{m}$  voxel size, 50 kV, no filter, 4,000 x 4,000 x 1,700 px – 10 mm sample size

## Pharmaceuticals & Packaging

- Measure coating thickness and distribution of active ingredients
- Check mechanical properties and defects
- Investigate pharmaceutical packaging up to a size of 7 cm x 7 cm x 7 cm
- Monitor and control the quality of metal and plastic components



3D rendering of paracetamol tablet. The coating has been color-coded for thickness  
5  $\mu\text{m}$  voxel size, 50 kV, no filter, 2,000 x 2,000 x 8,000 px – 11 mm sample size

# Available Now: 3D X-ray Vision for Everyone – X-ray Microtomography



X-ray micro-computed tomography ( $\mu$ CT) is one of the most advanced methods for getting 3D insights into samples of any material, any shape, and any size with little to no sample preparation.

Bruker microCT, a pioneer of  $\mu$ CT, has now made this technology easier and more accessible

for everyone by offering unparalleled 3D X-ray microscopy, all in the small size, Plug'n Analyze™ SKYSCAN 1272 high-resolution desktop  $\mu$ CT.

A single scan is all you need to reveal the complete internal 3D structure of your sample non-destructively.

# SKYSCAN 1272

## An Ingenious System & Powerful Desktop Solution for X-ray Microtomography

### **Genius-Mode™ or keep everything under user's control**

No prior  $\mu$ CT experience? No problem. The SKYSCAN 1272 offers automatic selection of parameters with Genius-Mode: magnification, energy, filter, and exposure time can all be optimized automatically with a single click. The intuitive interface gives more experienced users the freedom to adjust settings as desired for the ultimate scan.

### **Small footprint for the largest pictures: over 200 Megapixels in every slice**

Using a 16-megapixel X-ray camera in triple offset mode combined with the exclusively provided world's fastest reconstruction algorithm InstaRecon®, slices up to 14,456 x 14,456 pixels can be acquired. The large-area camera delivers high resolution over large scanning volumes, while keeping the working distance short for the best data quality. This huge image format allows zooming into any part of the 3D volume without rescanning the sample.

### **No hidden costs: a maintenance-free desktop $\mu$ CT scanner**

Our sealed X-ray tube allows running 24/7 without the frequent downtime required for changing a broken filament, saving you a lot of time and money.


### **One-stop shop: full in-house knowledge and global support**

Rely on Bruker's 20 years of experience in X-ray computed tomography and global team of experts for completely integrated hardware and software solutions. Come join our growing user family!






Various sample holders and in-situ stages



High resolution and ultimate ease-of-use for both the expert and novice



Smart solutions and design, such as the integrated vibration isolation, form a perfect whole



SKYSCAN 1272 high-resolution X-ray microtomograph



## Scanning 16 samples made easy!

Video at:  
[www.bruker.com/SKYSCAN1272-Video1](http://www.bruker.com/SKYSCAN1272-Video1)



Easy installation – just mount the sample changer on top of the scanner




Change samples at any time without interrupting an ongoing scan




Autodetection of new samples and status LEDs for every scan: ready, running, done



Rely on Genius-Mode for whatever task you may have!




Sample changer window



Scanning protocol window

- 1 Status display of all 16 positions
- 2 Automatic or user-selected parameters
- 3 All types of samples in the same tray

# Ready to Run 24/7 – as Automated as You Choose – Thanks to Genius-Mode™



SKYSCAN 1272 with automatic 16-position sample changer

## **SKYSCAN 1272 with 16-position sample changer and Genius-Mode means: multiple sample types, objects, and materials can be inter-mixed and scanned with varying protocols, automatically.**

The SKYSCAN 1272 with sample changer can be operated in three ways:

### **1) Fully automatic**

Simply load the sample changer, select “Auto” protocol, and then let the SKYSCAN 1272 take care of the rest! All scan settings are defined using Genius-Mode. Feel confident that your work is being done – all day, all night,

or over the weekend – with system-generated reports emailed directly to your inbox, including a link to access data remotely.

### **2) User selected**

Want more control? Individually adjust scan parameters for one, some, or all sixteen samples. Once all “Manual” protocols are defined, simply press “Start” to initiate the full batch.

### **3) Prior selection**

Streamline workflow by using the “Previous” command to assign the last settings.

Stay in charge, always. Because the sample changer operates outside the fully shielded X-ray chamber, a user can easily place a priority sample at the next position while another scan is still running.

# As Easy as 1-2-3: Let the Genius Work for You – Fully Automated

1.  
Moving to  
the **Best-Scan-Geometry™**


Thanks to SKYSCAN 1272's movable camera and its extra large X-ray beam opening, Genius-Mode finds the Best-Scan-Geometry – as compact as possible with the largest magnification – automatically.



Sample is too far away from the X-ray source, low magnification



Sample is too close and does not fit the field-of-view



Best-Scan-Geometry means: maximum magnification and minimum scan time due to the most compact setup


2.  
Finding the  
**Best-X-ray-Energy-Window™**

To find the perfect X-ray energy window, SKYSCAN 1272 automatically checks which of the six filters and X-ray energy best fits the sample's density in order to achieve the optimal image contrast.


Low Attenuation



50 kV, no filter



60 kV, Al 0.25 mm



70 kV, Al 0.5 mm



80 kV, Al 1 mm



90 kV, Al 0.5 + Cu 0.038 mm



100 kV, Cu 0.25 mm

High Attenuation



50 kV, no filter



60 kV, Al 0.25 mm



70 kV, Al 0.5 mm



80 kV, Al 1 mm




90 kV, Al 0.5 + Cu 0.038 mm



100 kV, Cu 0.25 mm

3.  
Starting with  
the **Best-Scan-Conditions™**

The SKYSCAN 1272 operating in Genius-Mode selects the best exposure time and rotation step automatically.



For highest resolution rely on the integrated micro-positioning stage and insert the sample manually



# SKYSCAN 1272 with Best-Scan-Geometry™ and Genius-Mode™

# Conventional Systems with Fixed Camera Position

## More Intensity and Speed

Best Position, Best Intensity, Best Speed



Fixed Position, Less Intensity, Less Speed

Yes




Variable Camera Position

No




Speed Factor



Moving both the sample and the large-format 16-MP camera as close as possible to the source increases the measured intensity dramatically. That's why SKYSCAN 1272 scans up to five times faster than conventional systems.

Conventional fixed camera systems are stuck at one scanning speed regardless of the sample size or magnification.

200+ MP  
Slice




Pharmaceutical tablet scanned in triple offset  
14,456 x 14,456 px with 0.8  $\mu$ m resolution, 11 mm sample size

# SKYSCAN 1272

## Over 200 Megapixels in Every Slice – Bigger is Better

Live optical camera



### Extremely large image format for more details and faster results


The SKYSCAN 1272 uses a unique family of specially developed, extremely large-format X-ray cameras, with fast readout and precise temperature stabilization for an optimal signal-to-noise ratio.

The camera can be moved horizontally to acquire projection images in three different offset positions – left, center, and right – which are automatically stitched together into a 14,456 x 3,280 pixel image, similar to a panoramic picture. The resulting reconstructed cross sections yield more than 200 million pixels in a single image!

### Sealed X-ray tube optimized for big samples and large-format imaging


By using a maintenance-free, sealed X-ray tube with a wide beam opening, the sample can be moved closer to the source while still maintaining a large field-of-view. Maximum X-ray flux is always available thanks to constant power output (10 W), independent of the filter or voltage setting, for maximum scanning speed.

The integrated regulation of the spot size, dependent on the settings, ensures a long lifetime and maintenance-free operation.




Analysis of the coating, 1,200 x 1,000 px

Analysis of distribution of active ingredients and defects in the matrix, 1,200 x 1,000 px







4

### X-ray Image

- Live view of X-ray projection images
- Easy switching between raw and background-corrected X-ray images
- Direct dimensional measurements

5

### Video Image

- Sample inspection with live optical camera
- Allows positioning of the sample for the highest resolution

6

### Control Bar

- Slider control for magnification, object position, and rotation
- Quick control of X-ray camera modes
- Easy filter selection to match sample attenuation



### Menu

- Simple, uncluttered menu for scanner control
- HELP database for additional information about features and functions

1

### Toolbar

- Natural left-to-right scanning workflow using clearly labeled icons
- Quick links to entire SKYSCAN software suite

2

### Scout View

- Full-length overview of entire sample
- Automatic stitching of oversized images
- Just click and drag at different vertical positions for batch scanning

3

### Just start your analysis!

Video at:  
[www.bruker.com/SKYSCAN1272-Video2](http://www.bruker.com/SKYSCAN1272-Video2)



## SKYSCAN 1272 Software Suite means Ease-of-Use plus Enjoy-your-Work

Intuitive, simple, yet powerful – the SKYSCAN 1272 control software is designed to inspire finding out what's inside. The whole screen, including all menus and icons, is laid out in a straightforward, left-to-right manner that even a first-time user will find intuitive enough to start scanning right away. All major functions can be performed with a single click, allowing researchers to focus on analyzing their samples, rather than finding buttons or navigating nested menus.

Overview of the SKYSCAN 1272 control software window  
Aluminum foam  
6  $\mu$ m voxel size, 90 kV, Al 0.5 + Cu 0.038 mm filter, 2,452 x 1,640 px

## Ultrafast reconstruction of large datasets


3D reconstruction from acquired 2D projection images used to require intensive computational power, and used to be the bottleneck of  $\mu$ CT analysis. The two approaches to overcome this barrier are: invest in more hardware – or use clever software.

SKYSCAN 1272 is supplied with InstaRecon® – the fastest reconstruction software available – exclusively offered by Bruker microCT.

InstaRecon® provides results up to 100 times faster than conventional reconstruction engines, operating on a single computer.


# Three, Two, One... Launch: The World's Fastest Software for Large-Format Reconstruction – InstaRecon®

Only with  
**InstaRecon®**  
can large-format images  
be handled  
in reasonable  
time



Cross section format [Pixel]	Cross sections in reconstructed volume [#]	Projections used for reconstruction [#]	1 PC & InstaRecon®	1 PC & 8 x GPU cluster acceleration	1 PC & GPU acceleration	1 PC
14,456 × 14,456	2,610	8,100	10 hours	19 hours	140 hours	1,440 hours*
8,000 × 8,000	2,495	2,157	43 min	68 min	275 min	3,552 min
4,000 × 4,000	2,255	1,990	9 min	19 min	62 min	913 min
2,000 × 2,000	1,229	996	80 sec	98 sec	247 sec	3,480 sec
1,000 × 1,000	615	499	12 sec	14 sec	21 sec	238 sec


\* time estimated



3 orthogonal slices through a wood sample in DATAVIEWER

## **DATAVIEWER** Slice-by-slice inspection of 3D volumes and 2D/3D image registration


DATAVIEWER allows inspection of the reconstructed volume using orthogonal slices in any direction. Objects can be rotated, repositioned, and resliced using their new orientation for more convenient visualization and saving of more efficient subvolumes. The software includes intuitive tools for measurement of 3D distances. 2D and 3D image registration enables the exact alignment of multiple scans of the same sample, acquired at different times.



Volume rendered wood sample, showing a color-coded pore size distribution by means of CTVOX

## **CTVOX** Realistic visualization by volume rendering


CTVOX is an easy-to-use volume rendering package that provides precise control of visualization parameters, ensuring a realistic representation of all types of samples. CTVOX offers intuitive manipulation of the point-of-view, virtual slicing through objects, and full control of light, shadow, and surface properties. Creating attractive cover images and movies that impress has never been so easy.



Analysis of local pore size in a wood sample by means of CTAN

## **CTAN** 2D/3D image analysis and processing

Built over two decades from direct feedback from scientists all over the world, CTAN is one of the most used programs for quantitative image analysis. This package includes an extensive number of tools for region-of-interest selection, image segmentation and 3D measurements. Using the comprehensive library of embedded plugins or user-customized protocols, quantifying complex microstructures such as porosity, thickness, orientation, and many other properties is easy. Simplify large study sets by batch analysis.



Surface rendering of a wood sample by means of CTVOL

## **CTVOL** Built-in surface rendering

Surface models can be visualized in CTVOL, a flexible 3D viewing environment. Volumes can be exported in STL format, to allow 3D printing of the acquired scan data or further use in CAD and modelling programs.

# Best Components, Superior Technology and Utmost Quality for Saving Energy, Time & Money



Maintenance-Free

**No Hidden Costs,**  
No Compressor,  
No Filaments

~99% Uptime

**Runs 24/7 for Years**  
and Will Never  
Let You Down

Low Power Consumption

**Saves 21,000 kWh**  
Electrical Energy  
per Year\*

The SKYSCAN 1272 provides top performance with peace of mind for years and will never let you down. No hidden costs of ownership, because our X-rays are "green" and the system is future-proof. Save maintenance, energy, and time – and in the end a lot of money!

\* 24 hours, 300 days

## SKYSCAN 1272 with Sealed X-ray Tube

System

Yes



Mainten



Maintenance-free

Power

Yes



Low CO<sub>2</sub>



90 W



## Systems with Open X-ray Tube

## SKYSCAN 1272 with Sealed X-ray Tube

## Systems with Open X-ray Tube




# Rely on the Only One-Stop Shop for $\mu$ CT and Become Part of the Bruker User Family



Bruker microCT Academy for education and training



CTVOX mobile app with full functionality



SKYSCAN 1272 desktop high-resolution  $\mu$ CT

## Hardware

- Systems, sample stages, computers and monitors from one supplier
- Fully calibrated and extensively tested hardware
- Direct installation and support from certified service engineers

## Software

- World's fastest reconstruction algorithm, InstaRecon<sup>®</sup>
- Powerful 3D analysis software and realistic 3D visualization
- Dedicated mobile app with full functionality and performance
- Multiple file formats for reporting and presentation
- Fully in-house developed software

## Experts

- Direct customer support and dedicated in-house experts
- Full system and software training
- Scientific support for applications and analytical tasks
- Newsletter with method training notes

**Get your CTVOX App & check out some samples!**

iOS

Android





Bruker microCT employs a team of researchers, engineers and technicians to provide cutting-edge desktop and laboratory  $\mu$ CT systems. From hardware to software, all of our experts work closely together and with customers to provide the best solution. Welcome to the only one-stop shop for X-ray microtomography.

By relying on a SKYSCAN 1272 you become part of the Bruker user family and benefit from the exchange of knowledge and experience. Bruker organizes annual  $\mu$ CT user meetings.

**We look forward to meeting you at our next get-together.**

**Get linked to the Bruker microCT Academy**


[www.bruker.com/products/microtomography/academy/academy.html](http://www.bruker.com/products/microtomography/academy/academy.html)



$\mu$ CT user meeting in Leuven, Belgium

## Overview of Features and Benefits

	Specification	Benefit
<b>X-ray source</b>	20 – 100 kV, 10W < 5 $\mu\text{m}$ spot size at 4 W	Covers a wide range of applications, from organics to metals
<b>Nominal resolution</b> (pixel size at maximum magnification)	< 0.35 $\mu\text{m}$ for 16 MP camera < 0.45 $\mu\text{m}$ for 11 MP camera	Detection of very small sample details
<b>X-ray camera</b>	16 MP, 4,904 x 3,280 px or 11 MP, 4,032 x 2,688 px 14-bit, cooled CCD fiber-optically coupled to scintillator	High resolution and large field-of-view
<b>Reconstructed volume</b> (after single scan)	up to 14,456 x 14,456 x 2,630 px for 16 MP camera up to 11,840 x 11,840 x 2,150 px for 11 MP camera	
<b>Sample size</b>	Max. $\varnothing$ 75 mm Max. height 70 mm	Fits small- to medium-sized objects
<b>Radiation safety</b>	< 1 $\mu\text{Sv/h}$ at any point on the instrument surface	Meets international safety requirements Easy installation
<b>Power supply</b>	100 – 240 V / 50 – 60 Hz	Standard wall socket Plug'n Analyze™
<b>Dimensions</b> (W x D x H)	116 cm x 52 cm x 33 cm 116 cm x 52 cm x 44 cm, with sample changer	Fits through standard doors Easy installation



**Bruker microCT**  
info.bmct@bruker.com

**Worldwide offices**  
bruker.com/baxs-offices

**Online information**  
bruker.com/microct

